


前额假体							单位 (CM)
		型号	A	B	C	D	用途
		K8010VS-1	3.5	6	0.3	3	隆前额手术
		K8010VS-2	4	6.5	0.3	3.3	隆前额手术
		K8010VS-3	4.5	7	0.3	3.8	隆前额手术
		K8010VS-4	5	8	0.35	3.8	隆前额手术
		K8010VS-5	5.5	9	0.35	4.7	隆前额手术
		K8010VS-6	6	10	0.35	5.2	隆前额手术
		K8010VS-7	6.5	11	0.4	5.7	隆前额手术
		K8010VS-8	7	12	0.4	6	隆前额手术
		K8010VS-9	7	13	0.4	6	隆前额手术
		K8010VS-10	7	14	0.5	6	隆前额手术